



PRODUCT SPECIFICATION

CATALOG NUMBER: G-107-C

| | |
|---------------|---|
| PRODUCT: | Goat anti-human IgE, affinity pure |
| FORM: | Liquid |
| PURITY: | Affinity purified |
| BUFFER: | 10 mM sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2 (PBS) |
| PRESERVATIVE: | 0.1% (w/v) sodium azide |
| STORAGE: | 2-8 degrees C |
| SPECIFICITY: | This antibody reacts with heavy chains on IgE, but not with other immunoglobulins in human serum. |