



PRODUCT SPECIFICATION

CATALOG NUMBER: G-112-C

| | |
|----------------------|--|
| PRODUCT DESCRIPTION: | Affinity Purified Goat anti-C reactive protein antibody |
| FORM: | Liquid |
| BUFFER: | 10 mM sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2 (PBS) |
| PRESERVATIVES: | 0.05% NaN ₃ |
| STORAGE TEMPERATURE: | -20°C long term; 2-8°C short term |
| HAZARD/BIOHAZARD: | This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product. |