



PRODUCT SPECIFICATION

CATALOG NUMBER: G13020

PRODUCT DESCRIPTION:	Goat anti-Human hCG Polyclonal, Beta, Affinity Purified
HOST ANIMAL:	Goat
IMMUNOGEN:	Human hCG Beta Subunit
FORM:	Liquid
BUFFER:	PBS, pH 7.4
PRESERVATIVE	0.1 % NaN ₃
STORAGE TEMPERATURE:	2-8°C short term; -20°C long term
ORIGIN:	U.S. Origin
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.