



PRODUCT SPECIFICATION

CATALOG NUMBER: G14020

PRODUCT DESCRIPTION:	Goat anti-Human HCG Polyclonal, Alpha, Affinity Purified
HOST ANIMAL:	Goat
IMMUNOGEN:	Human HCG Alpha
FORM:	Liquid
BUFFER:	0.02M PBS, buffer containing 0.15 M NaCl, pH 7.2
PRESERVATIVES:	0.1% NaN ₃
STORAGE TEMPERATURE:	2-8°C
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.