



PRODUCT SPECIFICATION

CATALOG NUMBER: G52230

PRODUCT DESCRIPTION:	Sheep anti Benzodiazepine Affinity Purified Polyclonal
IMMUNOGEN:	Benzodiazepine-BSA Conjugate
FORM:	Purified IgG from sheep serum
PURITY:	>90% by SDS-PAGE
PURIFICATION:	DEAE purified
BUFFER:	10 mM phosphate, 0.15M NaCl, pH 7.0
PRESERVATIVE:	0.02% sodium azide
STORAGE TEMPERATURE:	2-8°C short term; -20°C long term
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precaution. There is no known biohazard associated with this product.