



PRODUCT SPECIFICATION

CATALOG NUMBER: G82120

PRODUCT DESCRIPTION:	Goat anti-Human IgA Polyclonal, Alpha, Affinity Purified
FORM:	Liquid
PURITY:	The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
BUFFER:	0.01M Phosphate, 0.25M NaCl, pH 7.6
PRESERVATIVES:	None
STORAGE TEMPERATURE:	2-8°C under sterile conditions
SPECIFICITY:	Antibody reacts with heavy chain portion of IgA. This antibody does not react with heavy chains of human IgG or IgM. No antibody was detected against non-immunoglobulin serum proteins. The antibody may cross-react with IgA from other species.
COUNTRY OF ORIGIN:	USA
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.