



PRODUCT SPECIFICATION

CATALOG NUMBER: D-401-C-ABS

| | |
|----------------------|--|
| PRODUCT: | Donkey Anti-goat IgG (H & L) Antibody, Affinity Purified |
| FORM: | Liquid |
| PURITY: | Affinity purified |
| BUFFER: | 10 mM sodium phosphate, 0.85% (w/v) sodium chloride, pH 7.2 (PBS) |
| PRESERVATIVE: | 0.05% (w/v) sodium azide |
| STORAGE TEMPERATURE: | 2-8° C |
| SPECIFICITY: | This antibody reacts with goat IgG. No reactivity is observed to human IgG or rat IgG. |
| HAZARD/BIOHAZARD: | This material contains Sodium Azide as a preservative. Sodium azide may react with lead and copper plumbing to form explosive metal azides. Flush with large amounts of water during disposal. |