



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J13000**

PRODUCT DESCRIPTION:	HCG Beta Antigen, affinity pure
SOURCE:	Human pregnancy urine
FORM:	Lyophilized powder, lyophilized from 50 mM ammonium bicarbonate.
PURITY:	>98% by SDS PAGE
STORAGE TEMPERATURE:	2-8°C short term, -20°C long term
LEVEL OF CONTAMINANTS:	Whole HCG: 1.41% by Monobind ELISA FSH: 0.000006% by Bayer ACS 180 TSH: 0.0000023% by Bayer ACS 180 LH: 0.0000001% by Bayer ACS 180
HAZARD/BIOHAZARD:	Serum from donors tested negative for HBsAg, HBeAg, anti-HCV, anti-HIV 1&2, and HIVAg by FDA approved test methods. There is no known hazard associated with this product.