



**PRODUCT SPECIFICATION**

**CATALOG NUMBER: J13100**

PRODUCT DESCRIPTION:	HCG Antigen, High Purity	
SOURCE:	Human pregnancy urine	
FORM:	Lyophilized powder, lyophilized from 50 mM ammonium bicarbonate	
PURITY:	>95% by SDS PAGE	
IMMUNOLOGICAL ACTIVITY:	10,000 IU/MG	
LEVEL OF CONTAMINANTS:	hFSH:	0.0135%
	hLH	0.00381%
	hTSH:	0.00018%
STORAGE TEMPERATURE:	2-8°C short term -20°C long term	
HAZARD/BIOHAZARD:	Raw materials have been tested and found negative for HBsAg, anti-HCV antibodies, HAV antibodies, anti-HIV antibodies, and HIV antigen by FDA approved methods. There is no known hazard associated with this product.	