



PRODUCT SPECIFICATION

CATALOG NUMBER: J18010

PRODUCT DESCRIPTION:	FSH Antigen, standard grade
SOURCE:	Human Pituitary Glands
FORM:	Lyophilized powder, lyophilized from 50mM ammonium bicarbonate
PURITY:	>50% by SDS-PAGE
ACTIVITY:	1600 IU/mg
REFERENCE ASSAY:	FSH DPC Immulite 1000
REFERENCE PREP.:	WHO 2 nd IS 78/549
CONTAMINANTS:	LH: 0.0006% by Bayer ACS180 TSH: 0.000181% by Bayer ACS180
STORAGE TEMPERATURE:	-20° C long term 2-8° C short term
HAZARD/BIOHAZARD:	Donor material tested negative for HIV 1+2, HCV, HbsAg and Syphilis. There is no known hazard associated with this product.