



PRODUCT SPECIFICATION

CATALOG NUMBER: J18020

PRODUCT DESCRIPTION:	FSH Antigen, High Purity		
SOURCE:	Human Pituitary Glands		
FORM:	Lyophilized powder, lyophilized from 50mM ammonium bicarbonate		
PURITY:	>95% by SDS-PAGE		
REFERENCE ASSAY:	FSH DPC Immulite 1000		
REFERENCE PREP:	2 nd IS 78/549, WHO		
CONTAMINANTS:		Contaminant %	Contaminant Units
	LH	0.0014%	0.1190 IU/mg
	TSH	0.00494%	0.00042 IU/mg
	hGH	0.0711%	0.0007 mg/mg
	hPRL	0.0039%	0.000039 mg/mg
STORAGE TEMPERATURE:	2-8 C short term -20°C long term		
SPECIFIC ACTIVITY:	3800 IU/mg		
BIOHAZARD:	Tested negative for anti-HIV 1&2, anti-HCV, and HBsAg		