



PRODUCT SPECIFICATION

CATALOG NUMBER: J18050

PRODUCT DESCRIPTION:	Follicle Stimulating Hormone, FSH, Recombinant Human Antigen
SOURCE:	<i>Chinese Hamster Ovary Cells.</i>
FORM:	Lyophilized powder; lyophilized from a 1 mg/ml solution in PBS. Reconstitute in sterile deionized water at 0.1 mg/ml to 1 mg/ml, add 0.1% BSA or HSA to improve stability.
PURITY:	>98% by SDS-PAGE
STORAGE TEMPERATURE:	2-8°C short term, -20°C long term. After reconstitution, store antigen at -20°C. Avoid repeated freeze thaw cycles.
HAZARD/BIOHAZARD:	There is no known hazard or biohazard associated with this product.