



PRODUCT SPECIFICATION

CATALOG NUMBER: J31320

PRODUCT DESCRIPTION:	CKMB Antigen, Recombinant
SOURCE:	E. coli
FORM:	Frozen Liquid
PURIFICATION METHOD:	Combination of ion exchange chromatography and affinity column.
PURITY:	92.8% by SDS-PAGE
MOLECULAR WEIGHT:	~ 43 kD by SDS-PAGE for each subunit
BUFFER:	50 mM Tris-HCl buffer, 5 mM beta mercaptoethanol, pH 8.0, 1 mM EDTA, 5 mM disodium succinate
STABILIZER:	Glycerol
STORAGE TEMPERATURE:	-80°C
HAZARD/BIOHAZARD:	There is no known hazard or biohazard associated with this product.