



PRODUCT SPECIFICATION

CATALOG NUMBER: J37030

PRODUCT DESCRIPTION:	Recombinant Human Glycogen Phosphorylase BB (GPBB)
FORM:	Liquid
PURITY:	>90% by SDS-PAGE. Produced in <u>E. coli</u> , 97kD protein.
BUFFER:	25 mM beta glycerophosphate, 1 mM EDTA, 1 mM beta mercaptoethanol, pH 7.0. Glycerol added to 50%.
STORAGE TEMPERATURE:	-20°C or lower.