



PRODUCT SPECIFICATION

CATALOG NUMBER: J38010

| | |
|----------------------|---|
| PRODUCT DESCRIPTION: | Human D-Dimer |
| SOURCE: | Human plasma |
| FORM: | Frozen liquid |
| PURITY: | >90% by SDS PAGE |
| BUFFER: | 0.05 NaH ₂ PO ₄ buffer, pH 7.5 containing 0.15M NaCl and plus protease inhibitor cocktail. |
| PRESERVATIVE: | 0.09% NaN ₃ |
| STORAGE: | -15°C |
| IDENTITY: | One band at approx. 200kDa corresponding to MW of D-Dimer (4-20% PAGE) |
| HAZARD/BIOHAZARD: | This material contains sodium azide as a preservative, take appropriate precautions. Plasma from donors tested negative for HBsAg, and for anti-HIV-1, anti-HIV-2, anti-HCV, and syphilis by FDA approved test methods. |