



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J43040**

PRODUCT DESCRIPTION:	Cytomegalovirus (CMV) Glycoprotein B antigen, recombinant aa 11-76
SOURCE:	<u>E. coli</u>
PURIFICATION METHOD:	GS-4B Sepharose-affinity purification
PURITY:	>95% pure by 10% PAGE (coomassie staining)
BUFFER:	25 mM Tris-HCl, 1mM EDTA, 50% glycerol, pH 8.0
STORAGE TEMPERATURE:	Shipped at 2-8 °C Store at -20°C
SPECIFICITY:	Immunoreactive with sera of CMV-infected individuals.
STABILITY:	Five years frozen
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.