



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J43060**

PRODUCT DESCRIPTION:	CMV Pp28, Recombinant Antigen
SOURCE:	<u>E. coli</u>
PURIFICATION METHOD:	Sepharose-Derived
PURITY:	>95% pure by 10% PAGE (Coomassie staining)
BUFFER:	25 mM Tris-Hcl, 1mM EDTA, 50% glycerol, pH 7.2
STORAGE TEMPERATURE:	Shipped at ambient temperature. Upon arrival store at -20°C.
STABILITY:	Five years at -20°C; one month at room temperature.
SPECIFICITY:	Immunoreactive with sera of CMV-infected individuals.
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.