



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J43070**

|                      |   |
|----------------------|---|
| PRODUCT DESCRIPTION: | Cytomegalovirus (CMV) Pp38 recombinant antigen aa 117-373 |
| SOURCE:              | <u>E. coli</u>  |
| PURIFICATION METHOD: | GS-4B Sepharose-affinity purification                     |
| PURITY:              | >95% pure by SDS-PAGE (Coomassie staining)                |
| BUFFER:              | 25 mM Tris-HCl, 1mM EDTA, 50% glycerol, pH 7.2            |
| STORAGE TEMPERATURE: | 2-8 °C short term,<br>-20°C long term                     |
| STABILITY:           | Five years frozen   |
| SPECIFICITY:         | Immunoreactive with sera of CMV-infected individuals.     |
| HAZARD/BIOHAZARD:    | There is no known biohazard associated with this product. |