



PRODUCT SPECIFICATION

CATALOG NUMBER: J43080

PRODUCT DESCRIPTION:	Cytomegalovirus (CMV) Pp52 recombinant antigen aa 202-434
SOURCE:	<u>E. coli</u>
PURIFICATION METHOD:	Sepharose-derived purification
PURITY:	>95% pure by SDS-PAGE (Coomassie staining)
BUFFER:	50 mM Tris-HCl, 60 mM NaCl, 10 mM glutation,50% glycerol
STORAGE TEMPERATURE:	2-8 °C short term, -20°C long term
STABILITY:	Five years frozen
SPECIFICITY:	Immunoreactive with sera of CMV-infected individuals.
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.