



PRODUCT SPECIFICATION

CATALOG NUMBER: J43090

PRODUCT DESCRIPTION:	Cytomegalovirus (CMV) Pp65 antigen recombinant aa 297-510
SOURCE:	<u>E. coli</u>
PURIFICATION METHOD:	GS-4B sepharose-affinity purification
PURITY:	>95% pure by SDS-PAGE (Coomassie staining)
BUFFER:	25 mM Tris-HCl, 1 mM EDTA, 50% glycerol, pH 7.2
STORAGE TEMPERATURE:	2-8°C short term -20°C long term
SPECIFICITY:	Immunoreactive with sera of CMV-infected individuals.
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.