



PRODUCT SPECIFICATION

CATALOG NUMBER: J43510

PRODUCT DESCRIPTION:	EBV-VCA. Epstein Barr Virus, Viral Capsid Antigen
FORM:	Liquid
BUFFER:	20 mM Tris-HCl containing 3M Magnesium Chloride, pH 7.4
PRESERVATIVE:	None
STORAGE TEMPERATURE:	-70°C
PREPARATION:	Chronically infected cells are harvested by centrifugation and disrupted by sonication and cell debris is removed by centrifugation. This lysate is purified on an immunoaffinity column using a monoclonal antibody against gp125.
REACTIVITY:	This antigen preparation shows strong reactivity on Western Blots with a monoclonal antibody known to be specific for EBV gp 125.
BIOHAZARD:	All materials should be handled as if potentially infectious. Generally accepted laboratory practices appropriate for infectious materials should be employed when handling this product. There is no known biohazard associated with this product.