



PRODUCT SPECIFICATION

CATALOG NUMBER: J43540

PRODUCT DESCRIPTION:	EBV Mosaic EBNA1, Recombinant Antigen
SOURCE:	<u>E. coli</u>
PURIFICATION METHOD:	Sepharose-derived purification
PURITY:	>95% pure by 10% PAGE (Coomassie staining)
BUFFER:	50 mM Tris-HCl, 10mM glutation, 60 mM NaCl, 0.5% sarcosyl
STORAGE TEMPERATURE:	Shipped at ambient temperature. Upon arrival, store at -20°C.
STABILITY:	Five years frozen. One month in solution at room temperature.
SPECIFICITY:	Immunoreactive with sera of EBV-infected individuals.
HAZARD/BIOHAZARD:	There is no known hazard or biohazard associated with this product.