



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J43560**

PRODUCT DESCRIPTION:	EBV p23 Recombinant Antigen
SOURCE:	<u>E. coli</u>
PURIFICATION METHOD:	Sepharose-derived purification
PURITY:	>95% pure by SDS-PAGE (Coomassie staining)
BUFFER:	25 mM glycine, 50% glycerol, pH 9.6
STORAGE TEMPERATURE:	-20°C
SPECIFICITY:	Immunoreactive with sera of EBV-infected individuals.
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.