



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J44300**

PRODUCT DESCRIPTION:	Hepatitis A Antigen
FORMULATION:	Prepared from the supernatant of infected Wannenstapel and dialysed against PBS.
BUFFER:	Phosphate buffered saline
STORAGE TEMPERATURE:	-20°C
INACTIVATION:	Inactivated by incubating with formaldehyde for 11 days.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. The reagent has been inactivated and should contain no infectious material. There is no known biohazard associated with this product.