



## PRODUCT SPECIFICATION

CATALOG NUMBER: J44360

PRODUCT DESCRIPTION:	HAV VP1, Recombinant Antigen
SOURCE:	<u>E. coli</u>
PURIFICATION METHOD:	Inclusion bodies
PURITY:	90% pure by SDS-PAGE (Coomassie staining)
BUFFER:	10 mM CBB, 0.1% SDS, 50% glycerol, pH 9.6
STORAGE TEMPERATURE:	2-8°C short term, -20°C long term
STABILITY:	Five years at -20°C; one month at room temperature.
SPECIFICITY:	Immunoreactive with sera of HAV-infected individuals.
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.