



PRODUCT SPECIFICATION

CATALOG NUMBER: J44400

PRODUCT:	Recombinant Hepatitis B Core Antigen (rHBcAg)
SOURCE:	<u>E. coli</u>
FORM:	Liquid
PURIFICATION:	DEAE-52 and Sephacryl S-300
PURITY:	>80% by SDS-PAGE
BUFFER:	0.01 M PBS, 50% Glycerol, pH 7.4
PRESERVATIVE:	0.03% sodium azide
STORAGE TEMPERATURE:	2-8°C
REACTIVITY:	Reacts strongly with human anti-HBcAg positive serum