



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J47620**

PRODUCT DESCRIPTION:	HEV ORF2/ORF3, Recombinant Antigen
SOURCE:	<u>E. coli</u>
PURIFICATION METHOD:	GS-4B Sepharose-Affinity Purification
PURITY:	>95% pure by SDS-PAGE (Coomassie staining)
BUFFER:	25 mM Tris-Hcl, 1 mM EDTA, 0.5M urea, 50% glycerol
STORAGE TEMPERATURE:	2-8°C for shipping, -20°C long term
SPECIFICITY:	Immunoreactive with sera of HEV-infected individuals.
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.