



PRODUCT SPECIFICATION

CATALOG NUMBER: J61010

PRODUCT DESCRIPTION:	Alpha Fetoprotein
SOURCE:	Human cord serum
FORM:	Liquid
PURITY:	≥60%
BUFFER:	PBS, pH 7.3
PRESERVATIVES:	0.05% NaN ₃
STORAGE TEMPERATURE:	4°C short term. -20°C long term.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. Starting material chord serum tested and found to be negative for HBsAg, HIV 1 and HIV 2 antibodies, and HCV antibody. There is no known biohazard associated with this product.