



PRODUCT SPECIFICATION

CATALOG NUMBER: J62000

PRODUCT DESCRIPTION:	Human Carcinoembryonic Antigen (CEA), High Purity
SOURCE:	Human Liver
PURITY:	>98% by SDS-PAGE and electrophoresis
PURIFICATION:	Purified by immunoaffinity column, sterile filtered
BUFFER:	0.1 M PBS, pH 7.3
PRESERVATIVES:	None
CONTAMINANTS:	Specifically depleted from mitogenic impurities
BIOBURDEN:	Less than 100 CFU/ml
STORAGE TEMPERATURE:	2-8°C for shipping; -20°C long term storage Avoid freeze thaw cycles.
BIOHAZARD:	Material has been tested at donor level and found to be negative for anti-HCV, anti-HIV-1&2 and HBsAg by FDA approved assays. There is no known hazard associated with this product.