



## PRODUCT SPECIFICATION

**CATALOG NUMBER: J64000**

PRODUCT DESCRIPTION:	Beta-2-Microglobulin antigen, high purity
SOURCE:	Human urine
FORM:	Liquid
PURITY:	>90%
BUFFER:	40Mm PBS, pH 7.2, 150mM NaCl, & 1mM EDTA sodium salt
PRESERVATIVE:	0.05% NaN <sub>3</sub>
TOTAL PROTEIN:	determined by Biuret Method
STORAGE TEMPERATURE:	2 - 8°C, minimum shelf life 1 year
BIOHAZARD:	Donors tested and found negative by: HbsAg EIA Ortho 3 <sup>rd</sup> Gen & HBV DNA (PCR) Anti-HCV by Ortho 3 <sup>rd</sup> Gen & HCV RNA (PCR) Anti HIV 1 & 2 by Ortho 3 <sup>rd</sup> Gen & HIV DNA (PCR)