



PRODUCT SPECIFICATION

CATALOG NUMBER: J64010

PRODUCT DESCRIPTION:	Beta-2-Microglobulin antigen, standard grade
SOURCE:	Human urine
FORM:	Liquid aseptically filtered
PURITY:	Partially pure ~60%
BUFFER:	20mM Phosphate Buffered Saline (pH 8.0-8.2) with 10% glycerol
PRESERVATIVE:	0.15% sodium azide
STORAGE TEMPERATURE:	2 - 8°C
ANALYSIS:	Western Blot; SDS-PAGE
BIOHAZARD:	Donor tested HIV 1&2 antibody, HCV antibody, and HBsAg negative by FDA approved method, Viral DNA/RNA for HBV, HIV1 and HCV.