



PRODUCT SPECIFICATION

CATALOG NUMBER: J66100

PRODUCT DESCRIPTION:	CA 19-9 Antigen
BUFFER:	Liquid in 50 m M Tris, 150 mM NaCl, pH 7.4
PRESERVATIVE:	15mM sodium azide
CROSS REACTIVITY:	CA 125 11.69% CA 15-3 <1%
STORAGE TEMPERATURE:	2-8°C
BIOHAZARD:	All donor units have been tested and found negative for HBsAg, HIV I/II and HCV antibodies.
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.