



PRODUCT SPECIFICATION

CATALOG NUMBER: J66120

PRODUCT DESCRIPTION:	CA 19-9 Antigen, high purity
SOURCE:	Human Fluids
BUFFER:	Supplied Liquid in PBS, pH 7.4
PRESERVATIVE:	0.1% NaN ₃
PROTEIN CONCENTRATION:	0.073 OD Units
TOTAL PROTEIN METHOD:	OD@280
STORAGE TEMPERATURE:	2-8°C short term, -20°C long term
KNOWN CONTAMINANTS:	Ca 15-3 – 441 U/ml (< 0.1%) Ca 125 –4000 U/ml (3.5%)
BIOHAZARD:	This product has been tested and found negative for HIV 1&2 antibody, HbsAg, HCV antibody, and HIV 1-Ag.
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.