



PRODUCT SPECIFICATION

CATALOG NUMBER: J66310

PRODUCT DESCRIPTION:	CA 15-3 Breast Cancer Antigen
BUFFER:	25mM Tris, pH 7.4, with 140mM NaCl
PRESERVATIVE:	0.15% sodium azide
TOTAL PROTEIN:	3.32 OD units A280nm
CROSS REACTIVITY:	CA 125 <1% CA 19-9 <1% CEA – 6.009 ng/ml AFP <1%
STORAGE TEMPERATURE:	2-8°C
VIRAL TESTING:	Donor tested negative for HBV DNA, NCV RNA and HIV DNA by PCR
BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions.