



PRODUCT SPECIFICATION

CATALOG NUMBER: J75000

PRODUCT DESCRIPTION:	<u>Helicobacter pylori</u> Antigen, Part Pure
FORM:	In vitro purified proteins from H. Pylori
BUFFER:	15 mM phosphate, 0.15 M sodium chloride, pH 7.0
PRESERVATIVE:	0.1% sodium azide
STORAGE CONDITION:	short term 2-8°C; long term -20°C
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.