



PRODUCT SPECIFICATION

CATALOG NUMBER: J81200

PRODUCT DESCRIPTION:	Human Complement C4
SOURCE:	Normal human serum
FORM:	Single homogenous batch supplied frozen
BUFFER:	0.5M sodium phosphate buffer, pH 7.0
PRESERVATIVE:	0.09% NaN ₃
STORAGE TEMPERATURE:	-15°C
BIOHAZARD:	This product has been tested at serum donor level by FDA approved methods and found negative for HIV I & II antibodies and Hepatitis B surface antigen and Hepatitis C antibodies.