



PRODUCT SPECIFICATION

CATALOG NUMBER: J81610

PRODUCT DESCRIPTION:	Human C-Reactive Protein
SOURCE:	Human serum / plasma (tested and found negative for HIV I & II antibodies, HBs antigen and HCV antibody.
PURITY:	>99%
BUFFER:	20mM Tris buffer, pH 8.0, 0.28M NaCl, and 5mM CaCl ₂
PRESERVATIVE:	0.09% NaN ₃
FILTER:	Sterile filtered through a 0.2uM membrane.
STORAGE TEMPERATURE:	2 - 8°C
ANALYSIS:	Polyacrylamide gradient gel electrophoresis and cellulose acetate electrophoresis show one band only. 2-D IEP vs antisera to whole normal human sera shows single arc corresponding to CRP>
BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. Donor starting material tested and found negative for HIV 1&2 antibody, HCV antibody, and HBsAg.