



PRODUCT SPECIFICATION

CATALOG NUMBER: J81820

PRODUCT DESCRIPTION:	Ferritin, Human spleen
PURITY:	>96%
BUFFER:	0.05M phosphate buffer, pH 7.5, 0.15M NaCl
PRESERVATIVES:	0.09% NaN ₃
STORAGE TEMPERATURE:	2-8°C
TESTING:	Tested negative HIV I&II antibodies, Hepatitis B surface antigen, and Hepatitis C antibodies.
ANALYSIS:	Polyacrylamide gradient gel and cellulose acetate electrophoresis show one major band. Total protein determined by the Lowry method.
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.