



PRODUCT SPECIFICATION

CATALOG NUMBER: J82310

PRODUCT DESCRIPTION:	Human IgE Antigen (polyclonal) High Purity
SOURCE:	IgE from human serum
PURITY:	>95% (SDS-PAGE, IEP)
BUFFER:	0.1 M Tris, 0.2 M NaCl, pH 7.5
PRESERVATIVE:	0.1% NaN ₃
STORAGE:	2-8°C short term, -20°C long term
TESTING:	All donor serum tested negative for HbsAg, anti-HIV1/2, and anti-HCV.