



PRODUCT SPECIFICATION

CATALOG NUMBER: P81000

PRODUCT DESCRIPTION:	Ferritin Free Human Serum, Lipid Stripped
PRESERVATIVES:	None
TOTAL PROTEIN:	6.1 gm/dl
FERRITIN CONTENT:	<1.0 ng/ml by Elecsys
BIOCHEMICAL ASSAY:	Albumin 4.1 gm/dL Globulin 2.0 gm/dL Cholesterol 1.0 mg/dL Triglycerides <4.0 mg/dL Total Bilirubin <0.1 mg/dL Uric Acid 0.0 mg/dL Glucose <2.0 mg/dL Creatinine 0.0 mg/ml Phosphorus 2.1 mg/dl Alk Phos 47.0 U/L Calcium <0.2 mg/ml Sodium 139.0 mEq/L Potassium 2.11 mEq/L Chloride 97.7 mEq/L ALT 6.0 U/L AST 14.0 U/L GGT 11.0 U/L
MICROBIOLOGY ASSAY:	CFU/mL <100 at bottling
STORAGE TEMPERATURE:	-20°C
BIOHAZARD:	Each human blood donor Unit of the above referenced product has been tested by FDA approved methods and found to be non reactive for HIV-1/2 Antibody, HCV antibody, HIV Ag, and HBsAg. The above referenced product was only formed with these donors.