



PRODUCT SPECIFICATION

CATALOG NUMBER: T82323

PRODUCT DESCRIPTION:	Goat anti-human IgE, affinity purified-HRP conjugate
BUFFER:	PBS containing 0.2% BSA and 0.1% Pro-Clean 400
MOLAR ENZYME/ANTIBODY PROTEIN RATIO:	4:1
DILUTION FACTOR:	1:10,000-1:100,000 (ELISA)
STORAGE:	2-8°C
COUNTRY OF ORIGIN:	USA
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.