



PRODUCT SPECIFICATION

CATALOG NUMBER: T82422

PRODUCT DESCRIPTION:	Affinity Purified Goat Anti Human IgG-Fc Alkaline Phosphatase Conjugate
FORM:	Liquid
BUFFER:	1.0MM MgCl ₂ , 0.1MM ZnCl ₂
STORAGE TEMPERATURE:	2-8°C
COUNTRY OF ORIGIN:	USA
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.