



## PRODUCT SPECIFICATION

**CATALOG NUMBER: T82740**

PRODUCT DESCRIPTION:	Affinity Pure Goat anti Human IgA, alpha chain specific, Alkaline Phosphatase conjugate
FORM:	Lyophilized powder
BUFFER:	Lyophilized in 0.01 M Tris-HCl, 0.25 M NaCl, pH 8.0, with 15 mg/ml BSA added as stabilizer.
PRESERVATIVE:	0.05% sodium azide
RECONSTITUTION:	Reconstitute with 1.0 ml distilled water, let stand at room temp for 2 hours and centrifuge if not clear.
SUGGESTED DILUTION:	1/5000-1/50,000 for ELISA and Western blots
STORAGE TEMPERATURE:	2-8°C until reconstitution. Store reconstituted antibody at 2-8°C short term. For long term storage add Glycerol to 50% and store at -20°C.
SPECIFICITY:	Reacts with human IgA but not human IgG or human IgM. May cross react with IgA from other species