



PRODUCT SPECIFICATION

CATALOG NUMBER: T83460

PRODUCT DESCRIPTION:	Donkey anti Goat IgG H&L Conjugate, FITC
FORM:	Freeze dried powder
BUFFER:	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6
PRESERVATIVE:	0.05% Sodium Azide
SUGGESTED DILUTION:	1:50 – 1:200 for most applications
FLUOROPHORE:	Fluorescein-5-isothiocyanate (FITC-isomer-1) Amax = 492 nm; Emax = 520 nm
FLUOROPHORE / PROTEIN:	8.6 ug/mg; 3.2 moles FITC per mole IgG
STABILIZER:	15 mg/mL Bovine Serum Albumin (IgG-Free, Protease-Free)
STORAGE:	2-8 ^o C short term, -20 ^o C long term
COUNTRY OF ORIGIN:	USA
HAZARD/BIOHAZARD:	This material contains sodium azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.
NOTE:	For in vitro use only, not for diagnostic or therapeutic use.