



PRODUCT SPECIFICATION

CATALOG NUMBER: T83822

PRODUCT DESCRIPTION:	Rabbit Anti-Goat IgG (H&L)-Alkaline Phosphatase Conjugated
FORM:	Liquid
PURIFICATION:	Isolated by affinity chromatography using antigen coupled to agarose beads and conjugated to alkaline phosphatase.
BUFFER:	50 mM HEPES, 0.1 M NaCl, 0.1 mM MgCl ₂ , 0.1 mM ZnCl ₂ , pH 7.1
PRESERVATIVE:	0.1% Sodium Azide
STABILIZING AGENT:	0.2% BSA
STORAGE:	2-8°C
HAZARD/BIOHAZARD:	This material contains Sodium Azide as a preservative, take appropriate precautions. There is no known biohazard associated with this product.