



## PRODUCT SPECIFICATION

**CATALOG NUMBER: V52000**

PRODUCT DESCRIPTION:	Benzoylcegonine-BSA Antigen Conjugate
FORM:	Salt free lyophilized protein
RECONSTITUTION:	Reconstitute in distilled water
METHOD LINK:	Activated Ester
STORAGE:	After reconstitution 2-8°C short term and -20°C long term. Avoid repeated freeze-thaw cycles.
HAZARD/BIOHAZARD:	There is no known biohazard associated with this product.